


8. Intraoral

8.1 Surgical Hooks

8.2 Lingual Accessories

8.1 Surgical Hooks

T-Bar Surgical Split 10/pk



Made by MIM technique and highest quality Stainless Steel, the high strength ball hooks easily crimp to archwires, in or out of the mouth. For use when extra traction is needed. Split ball hooks provide convenient positioning without removing archwires.


Our Ball Hook Crimpable pliers is to be suggested.



Slot	Part Code	length	Part Code	length
Square	F730-000	2.1mm	F731-000	5.97mm
Universal	F740-000	2.1mm	F741-000	5.97mm
Round	F750-000	2.1mm	F751-000	5.97mm

Slide-on Stops 10/pk

Slot	Part Code	length
Square	F730-100	2.1mm
Universal	F740-100	2.1mm
Round	F750-100	2.1mm



Split Ball Hook and Shepherd Hook 10/pk

Slot	Part Code	length		
Square	F730-101			
Shepherd Hook	8mm	10/pk	10mm	10/pk
Square	F750-101	Right	F750-201	Right
Square	F750-201	Left	F750-202	Left


8.2 Lingual Accessories

Palatal Arch 1/pk

Made of highly polished stainless steel, it offer ultimately durability and is compatible with all sterilization. Sold in set of four or individually.


NO.	Part Code	No	Part Code
37#	F600-100	47#	F600-600
39#	F600-200	49#	F600-700
41#	F600-300	51#	F600-800
43#	F600-400	53#	F600-900
45#	F600-500	55#	F600-1000

Expansion Screw 1/pk



Type	Guide Pins	Elongation	Length	Width	Thickness	Part No.
Standard	2	7mm	14mm	10mm	3.2mm	F870-000
Trestle	0	5mm/12mm	22mm	10mm	3.2mm	F880-000

Springs



	Clockwise	Counterclockwise
Rotation Spring	W320-800	W320-900
Uprighting Spring	W330-800	W33-900
Seperating Springs	W310-000	.05mm

Lingual Retainer



Upper		Lower	
38#	F210-100	29#	F220-100
40#	F210-200	31#	F220-200
42#	F210-300	33#	F220-300
45#	F210-400	35#	F220-400
47#	F210-500	37#	F220-500